

VANTAGE

Ordering Information



Product	Catalog Number	Description
VENTANA VANTAGE Workflow Solution	For further information contact your local sales representative	The VANTAGE workflow solution suite is designed to help improve productivity and monitor patient samples in anatomic pathology. This solution incorporates Lean Sigma principles and provides comprehensive specimen tracking, greater quality control and unparalleled informatics
Central Management	1339700	This networked software product provides single-point management of multiple advanced staining systems. Customers with multiple hosts can share reagents and protocols. Single-host customers can add a second host and start multiple runs at the same time.
VENTANA LIS Connectivity Interface	1339800	This networked product provides an interface between a VENTANA staining system and the laboratory information system (LIS).
VANTAGE Slide Barcode Label	1749300	Roll of 1700 blank slide labels. These slide labels have been validated for use on all Ventana staining platforms.
VANTAGE Zebra GX430t Printer Ribbon Roll	1747600	High-performance wax ribbon for use on the Zebra GX 430 label printers. This ribbon will produce approximately 2880 VANTAGE Slide Barcode labels.
VANTAGE Overlay	1749400	Roll of 3750 clear slide overlays. These overlays are required for use on the BenchMark Series platforms, when used with the VANTAGE Slide Barcode Label.